



In House Ocular Pack

AMS12947

CONTENTS:

2	BLADE SCALPEL SAFETY #15	1	COVER BACKTABLE 44 X 44 REIN	1	CUP MEDICINE 2 OZ CLEAR
1	DRAPE OPHTHALMIC BILATERAL	P 2	LAB VACUTAINER 6ML LAV K2	P 1	LAB VACUTAINER 10ML RED
P 2	LAB VACUTAINER 8.5ML W/CLOT	P 1	LAB VACUTAINER 6ML	2	GOWN SLEEVE NS
10	LABEL CUSTOM	5	LABEL CUSTOM	4	NEEDLE ECLIPSE 18 X 1.5
10	SPONGE 6 WOOD COTTON TIP	4	SPONGE COTTON BALL LG	10	SPONGE 4 X 4 X 12 PLAIN
1	SYRINGE 3CC L/L	4	SYRINGE 10CC L/L	P 1	TOWEL ABS 15 X 22 2/PK
1	TRAY PREP 9 X 5 X 2	2	WRAP CSR 24 X 24 CUSTOM	2	WRAP CSR 30 X 30 CUSTOM
5	FASTENER ZIP TIE 11.5"	1	BAG ZIPLOCK 10 X 12 2ML	1	BAG ZIPLOCK 12 X 15 2ML
1	BAG ZIPLOCK 4 X 6 2ML	5	BAG POLY 8 X 18 4MIL	P 2	BAG BIOHAZARD ZIPLOCK 8 X 10
2	PREP ALCOHOL PAD MD	2	PREP IODOPHOR PVP PAD	3	PREP PVP SWABSTICK 3'S
P 1	PREP SOLUTION EYE WASH 4 OZ	P 1	PREP SOLUTION OPHTHALMIC	P 1	ENVELOPE PACKING LIST
1	TAG DONOR 3 X 6 WHITE	P 2	CORNEAL VIEWING CHAMBER	1	PENLIGHT DISPOSABLE
P 2	EYE CAP LG OVAL FLESH COLOR	P 1	SALINE IRRIGATION 0.9% 100ML		

P - Item is Outside Pack



LOT:



ASSEMBLED IN THE U.S.A

Contents of the unopened, undamaged package are sterile. Do not re-sterilize. Contents of foil packets, bottles or ampules have not been terminally sterilized by AMS. Federal law (USA) restricts this device to sale by or on the order of a physician. Some components may be of Foreign Origin. Store at a controlled room temperature 15°C-30°C (59°F-86°F). Open using Hospital's Standard Practices. To prevent damage, it is recommended to store upright on a non-abrasive surface.