



STERILE EO

2D
BARCODE

Closure Kit AMS13043

CONTENTS:

- | | | |
|---------------------------|-------------------------------|-------------------------------|
| 1 BANDAGE EYE PAD | 1 BLADE #15 CS | 1 CAUTERY BLADE ELECTROLASE |
| 1 COVER LIGHT HANDLE FLEX | 1 COUNTER NEEDLE 10 CNT | 5 SPONGE 6 WOOD COTTON TIP |
| 3 SPONGE 8 OB POLY PROP | 1 DRESSING NON-ADHERENT 2 X 3 | 20 SPONGE NON WOVEN 4 X 4 X 4 |
| 3 TOWEL OR BLUE | 1 TRAY PREP 9 X 5 X 2 | 1 WRAP 28 X 28 WHITE TPT |
| 1 BAG ZIPLOCK 8 X 10 2ML | 1 SHEATH HANDPIECE | P 1 SALINE .9% 110ML CUP |

P - Item is Outside Pack



LOT:



LINEAR BARCODE

ASSEMBLED IN THE U.S.A

Contents of the unopened, undamaged package are sterile. Do not re-sterilize. Contents of foil packets, bottles or ampules have not been terminally sterilized by AMS. Federal law (USA) restricts this device to sale by or on the order of a physician. Some components may be of Foreign Origin. Store at a controlled room temperature 15°C-30°C (59°F-86°F). Open using Hospital's Standard Practices. To prevent damage, it is recommended to store upright on a non-abrasive surface.

ALIGNED MEDICAL SOLUTIONS | 1602 4TH AVE NORTH | BILLINGS MT 59101 | TEL 800-637-7056 | FAX 406-256-9875